

**SEKTION H**  
HNO



**SECTION H**  
ENT

**SECCIÓN H**  
ORL

<b>Inhaltsverzeichnis</b>	<b>Contents</b>	<b>Contenido</b>		<b>Seite, Page, Página</b>
<b>Sektion H</b>	<b>Section H</b>	<b>Sección H</b>		
Kehlkopfspiegel, Trachealdilatatoren, -häkchen, Trachealkanülen	Laryngeal Mirrors, Tracheal Dilators, Hooks and Tubes	Espejos laringeos, dilatadores, ganchos y cánulas traqueales	<b>HA</b>	H.01 - H.04
Zungendrucker, Mundsperrer, Zungenzange	Tongue Depressors, Mouth Gags, Tongue Forceps	Abaja-lenguas, Abre-bocas, Pinza tira-lenguas	<b>HC</b>	H.05 - H.07
Ligaturnadeln, Dissektoren, Tonsillenschnürer, -zangen, Conchotome, Ringmesser	Ligature Needles, Dissectors, Tonsil Snares, Tonsil Forceps, Punches, Adenoid Currettes	Agujas para ligaduras, disecadores, Cierra-nudos, Pinzas amígdales, Cuchillos adenoides	<b>HD</b>	H.08 - H.11
Nasenspekula	Nasal Specula	Espéculos nasales	<b>HG</b>	H.12 - H.13
Antrumküretten, Septum- elevatoren, Schwingmesser, Nasenpolypenschnürer, Septumhohlmeißel und -osteotome, Nasenraspeln, -sägen, Spühlkatheter	Antrum Currettes, Septum Elevators, Swivel Knives, Nasal Snares, Septum Gouges, Osteotomes, Nasal Rasps and Saws, Cannulas	Curetas para senos maxilares, elevadores y bisturíes para el septo, cierra-nudos para pólipos nasales, gubias y cinceles para el septo, raspas y sierras nasales, Canulas	<b>HH</b>	H.14 - H.20
Nasenpolypenzangen, Nasenconchotome, Kieferhöhlenstanzen, Septumzangen, -scheren	Nasal Polypus Forceps, Nasal Cutting Forceps, Antrum Punches, Septum Forceps and Scissors	Pinzas para pólipos nasales, pinzas cortantes nasales y para senos maxilares, pinzas y tijeras para el septo	<b>HK</b>	H.21 - H.24
Ohrtrichter, Stimmgabeln, Ohrkatheter, -spritzen, Ohrpolypenschnürer, -küretten, -schlingen, Parazentese- nadeln, Ohrhäkchen	Ear Specula, Tuning Forks, Eustachian Catheters, Ear Syringes, Ear Polypus Snares, Ear Currettes and Loops, Paracentesis Needles, Ear Hooks	Espéculos auriculares, diapaso- nes, tubos y jeringas auriculares, cierra-nudos auriculares, curetas y asas auriculares, agujas para para-centesis, ganchos auricula- res	<b>HO</b>	H.25 - H.31
Ohrpolypenzängchen, Fremdkörperzängchen, Ohrzängchen, Knochenstanzen	Ear Polypus Forceps, Foreign Body Forceps, Ear Forceps, Bone Punches	Pinzas para pólipos auriculares, Pinzas para cuerpos extraños, Pinzas auriculares, Punzonadoras de huesos	<b>HS</b>	H.32 - H.34
Setvorschläge	Set Suggestions	Proposición de Sets	<b>ZH</b>	H.35 - H.38

**Kehlkopfspiegel**  
**Laryngeal Mirrors**  
**Espejos laringeos**



mit Griff - with handle - con mango

<b>HA 101-08</b>	# 00	ø 8 mm
<b>HA 101-10</b>	# 0	ø 10 mm
<b>HA 101-12</b>	# 1	ø 12 mm
<b>HA 101-14</b>	# 2	ø 14 mm
<b>HA 101-16</b>	# 3	ø 16 mm
<b>HA 101-18</b>	# 4	ø 18 mm
<b>HA 101-20</b>	# 5	ø 20 mm
<b>HA 101-22</b>	# 6	ø 22 mm
<b>HA 101-24</b>	# 7	ø 24 mm
<b>HA 101-26</b>	# 8	ø 26 mm
<b>HA 101-28</b>	# 9	ø 28 mm
<b>HA 101-30</b>	# 10	ø 30 mm



ohne Griff - without handle - sin mango

<b>HA 105-08</b>	# 00	ø 8 mm
<b>HA 105-10</b>	# 0	ø 10 mm
<b>HA 105-12</b>	# 1	ø 12 mm
<b>HA 105-14</b>	# 2	ø 14 mm
<b>HA 105-16</b>	# 3	ø 16 mm
<b>HA 105-18</b>	# 4	ø 18 mm
<b>HA 105-20</b>	# 5	ø 20 mm
<b>HA 105-22</b>	# 6	ø 22 mm
<b>HA 105-24</b>	# 7	ø 24 mm
<b>HA 105-26</b>	# 8	ø 26 mm
<b>HA 105-28</b>	# 9	ø 28 mm
<b>HA 105-30</b>	# 10	ø 30 mm




Griff für Kehlkopfspiegel  
Handle for laryngeal mirrors  
Mango para espejos laríngeos  
**HA 110-10** 10 cm / 4"

**Trachealdilatatoren und -häkchen**  
**Tracheal Dilators and Hooks**  
**Dilatadores y ganchos traqueales**



**Laborde**  
**HA 121-14** 14 cm / 5½"



**Trousseau**  
**HA 151-14** 14 cm / 5½"




**Jackson**  
**HA 214-15** 15 cm / 6"



**Bose**  
 scharf - sharp - agudo

rechts gebogen, für die linke Hand  
 curved to right, for the left hand  
 curva a derecha, para la mano izquierda  
**HA 234-16** 16 cm / 6¼"



links gebogen, für die rechte Hand  
 curved to left, for the right hand  
 curva a izquierda, para la mano derecha  
**HA 235-16** 16 cm / 6¼"


**Trachealhäkchen**  
**Tracheal Hooks**  
**Ganchos traqueales**




**Iterson**

stumpf - blunt - romo

**HA 240-16** 16 cm / 6¼"



**HA 240-16**



**HA 244-16**

scharf - sharp - agudo

**HA 244-16** 16 cm / 6¼"



**Lukens**

**HA 251-17** 17 cm / 6¾"



**Jackson**

**HA 261-17** 17 cm / 6¾"

**Trachealkanülen**  
**Tracheal Tubes**  
**Cánulas traqueales**



**Luer**  
*Normalradius*  
*Normal radius*  
*Radio normal*

Neusilber - German silver - plata nueva

<b>HA 420-04</b>	# 0	ø 4 mm
<b>HA 420-05</b>	# 1	ø 5 mm
<b>HA 420-06</b>	# 2	ø 6 mm
<b>HA 420-07</b>	# 3	ø 7 mm
<b>HA 420-08</b>	# 4	ø 8 mm
<b>HA 420-09</b>	# 5	ø 9 mm
<b>HA 420-10</b>	# 6	ø 10 mm
<b>HA 420-11</b>	# 7	ø 11 mm
<b>HA 420-12</b>	# 8	ø 12 mm
<b>HA 420-13</b>	# 9	ø 13 mm



**Chevalier-Jackson**  
*Normalradius*  
*Normal radius*  
*Radio normal*

mit Obturator - with obturator - con obturador

Neusilber - German silver - plata nueva

<b>HA 425-04</b>	# 0	ø 4 mm
<b>HA 425-05</b>	# 1	ø 5 mm
<b>HA 425-06</b>	# 2	ø 6 mm
<b>HA 425-07</b>	# 3	ø 7 mm
<b>HA 425-08</b>	# 4	ø 8 mm
<b>HA 425-09</b>	# 5	ø 9 mm
<b>HA 425-10</b>	# 6	ø 10 mm
<b>HA 425-11</b>	# 7	ø 11 mm
<b>HA 425-12</b>	# 8	ø 12 mm
<b>HA 425-13</b>	# 9	ø 13 mm

**Zungendrucker**  
**Tongue Depressors**  
**Abaja-lenguas**




**Brünings**  
HC 070-19 19 cm / 7½"



**Wieder**  
HC 080-14 14 cm / 5½"  
HC 080-15 15 cm / 6"



**Hartmann**  
HC 083-15 15 cm / 6"




**Tobold**  
HC 101-14 14,5 cm / 5¾"

**Mundsperrer**  
**Mouth Gags**  
**Abre-bocas**



**Heister**  
HC 110-13 13 cm / 5"



**Molt**  
HC 121-14 14 cm / 5½"



**Denhart**  
HC 131-12 12 cm / 4¾"



**Jennings**  
HC 170-09 9 cm / 3½"  
HC 170-11 11 cm / 4¼"  
HC 170-13 13 cm / 5"  
HC 170-15 15 cm / 6"



**Mundsperrer, Zungenzange**  
**Mouth Gags, Tongue Forceps**  
**Abre-bocas, Pinza tira-lenguas**



HC 302-01



HC 302-02



HC 302-03



HC 302-04



HC 302-05



**Davis-Boyle**

komplett - complete - completo

HC 300-00

bestehend aus - including - compuesto de:

Sperrer allein - frame only - cuadro solo

HC 301-01

*Maße: Tiefe x Breite*

*Dimensions: Depth x Width*

*Dimensiones: Profundidad x Anchura*

Zungendrucker - tongue plates - abaja-lenguas

HC 302-01 26 mm x 21 mm

HC 302-02 65 mm x 23 mm

HC 302-03 80 mm x 24 mm

HC 302-04 93 mm x 25 mm

HC 302-05 108 mm x 26 mm



**Collin**

HC 700-17 17 cm / 6¾"

**Ligaturnadeln, Dissektoren, Tonsillenzangen**

**Ligature Needles, Dissectors, Tonsil Forceps**

**Agujas para ligaduras, disecadores, Pinzas amígdalas**



**Hurd**

links gebogen, für die rechte Hand  
curved to left, for the right hand  
curva a izquierda, para la mano derecha  
**HD 101-21** 21 cm / 8¼"

rechts gebogen, für die linke Hand  
curved to right, for the left hand  
curva a derecha, para la mano izquierda  
**HD 103-21** 21 cm / 8¼"



**Hurd**


**HD 171-22** 22,5 cm / 9"




**Henke**

**HD 191-23** 23 cm / 9"

**HD 193-23** 23 cm / 9"



**HD 191-23**



**HD 193-23**



**White**

gebogen - curved - curva  
**HD 431-18** 18 cm / 7"  
**HD 431-23** 23 cm / 9"

**Tonsillenschnürer**  
**Tonsil Snares**  
**Cierra-nudos amígdalas**



**Eves**

ohne Raste - without ratchet - sin cremallera  
**HD 470-28** 28 cm / 11"




**Eves**

mit Raste - with ratchet - con cremallera  
**HD 472-28** 28 cm / 11"



**Bruenings**

mit Raste - with ratchet - con cremallera  
**HD 490-28** 28 cm / 11"




Schnürdraht, geformt, Dtzd.

Snare wire, shaped, doz.

Alambre, formado, dc.

**HD 495-30** 30 mm x 100 mm



Schnürdrahtschlingen, Packung zu 100


Snare wire loops, package of 100

Lazos de alambre, paquete de 100


**HD 495-40** 40 mm

**HD 495-50** 50 mm

**Tonsillenschnürer und Conchotome**  
**Tonsil Snare and Punches**  
**Cierra-nudos y pinzas cortantes amígdalas**



**Tying**  
komplett, mit 1 geraden Aufsatz  
complete, with 1 straight tip  
completo, con 1 guía-asa recta  
**HD 500-01**




**Hartmann**  
14 cm / 5½"  
**HD 530-13**    ø 13 mm  
**HD 530-15**    ø 15 mm  
**HD 530-17**    ø 17 mm

**Ringmesser**  
**Adenoid Currettes**  
**Cuchillos adenoides**




**Beckmann**  
gerade - straight - recto  
22 cm / 8<sup>3</sup>/<sub>4</sub>"


- HD 551-01 # 1
- HD 551-02 # 2
- HD 551-03 # 3
- HD 551-04 # 4
- HD 551-05 # 5
- HD 551-06 # 6




HD 551-01




HD 551-02




HD 551-03



HD 551-04



HD 551-05




HD 551-06




**Barnhill**  
gerade - straight - recto  
22 cm / 8<sup>3</sup>/<sub>4</sub>"


- HD 561-01 # 0
- HD 561-02 # 1
- HD 561-03 # 2
- HD 561-04 # 3
- HD 561-05 # 4




HD 561-01




HD 561-02



HD 561-03



HD 561-04



HD 561-05

**Nasenspekulas**  
**Specula Nasal**  
**Espéculos nasales**



**Modell Wien**  
14,5 cm / 5¾"  
**HG 160-14** 30 mm



**Tieck-Halle**  
13,5 cm / 5¼"  
**HG 170-13** 17 mm



**Hartmann-Halle**  
15 cm / 6"  
**HG 200-01** 25 mm  
**HG 200-02** 30 mm  
**HG 200-03** 35 mm

**Nasenspekulas**  
**Specula Nasal**  
**Espéculos nasales**



**Killian**

13 cm / 5"

**HG 230-35** 35 mm

**HG 230-50** 50 mm

**HG 230-75** 75 mm

**HG 230-90** 90 mm



**Cottle**

13 cm / 5"

**HG 240-35** 35 mm

**HG 240-50** 50 mm

**HG 240-75** 75 mm

**HG 240-90** 90 mm



**Cottle**

schwarz verchromt  
black chrome-plated  
cromado negro

14 cm / 5½"

**HG 250-45** 45 mm

**HG 250-75** 75 mm

**Antrumküretten**

**Antrum Currettes**

**Curetas para senos maxilares**



**Pratt**

doppelendig - double-ended - dos extremos  
ø 5 mm

**HH 171-23** 23 cm / 9"



**Faulkner**

doppelendig - double-ended - dos extremos  
ø 5 mm / 8 mm

**HH 191-22** 22 cm / 8¾"



**Coakley**

17 cm / 6¾"

**HH 201-01** # 1

**HH 201-02** # 2

**HH 201-03** # 3

**HH 201-04** # 4

**HH 201-05** # 5

**HH 201-06** # 6



**HH 201-01**



**HH 201-02**



**HH 201-03**



**HH 201-04**



**HH 201-05**



**HH 201-06**



**Septumelevatorien, Schwingmesser**  
**Septum Elevators, Swivel Knives**  
**Elevadores y bisturías para el septo**



**Cottle**

doppelendig - double-ended - dos extremos

Breiten - widths - anchuras 4,5 mm / 4,5 mm  
**HH 230-22** 22,5 cm / 9"



**Cottle**

doppelendig - double-ended - dos extremos

Breiten - widths - anchuras 4,5 mm / 4,5 mm  
**HH 250-23** 23 cm / 9"



**Freer**

doppelendig, scharf / stumpf  
double-ended, sharp / blunt  
dos extremos, cortante / romo

Breiten - widths - anchuras 5 mm / 5 mm  
**HH 315-18** 18 cm / 7"




**Ballenger**

20 cm / 8"


**HH 400-03** 3 mm

**HH 400-04** 4 mm


**HH 400-05** 5 mm



**HH 400-03**



**HH 400-04**



**HH 400-05**

## Nasenpolypenschnürer

### Nasal Snares

### Cierra-nudos para pólipos nasales



Schnürdraht, Ring mit 10 m

Snare wire, coil of 10 m

Alambre, rollo de 10 m

**HH 405-03** ø 0,3 mm

**HH 405-04** ø 0,4 mm



Schnürdraht, geformt, Dtzd.

Snare wire, shaped, doz.

Alambre, formado, dc.

**HH 407-03** ø 0,3 mm



#### Krause

1 Aufsatz - 1 working end - 1 punta

**HH 411-26** 26 cm / 10¼"



#### Krause-Voss

2 Aufsätze - 2 working ends - 2 puntas

**HH 421-26** 26 cm / 10¼"




#### Lange-Wilde


1 Aufsatz - 1 working end - 1 punta

**HH 431-26** 26 cm / 10¼"

## Septumhohlmeißel und -osteotome Septum Gouges and Osteotomes Gubias y cinceles para el septo



HH 483-16



HH 485-16

### Kilian

Hohlmeißel - Gouges - Gubias  
bajonettförmig - bayonet shaped - forma bayoneta

16 cm / 6 1/4"

HH 483-16 5 mm

HH 485-16 5 mm



### Cottle

Osteotome - Osteotomes - Osteótomos  
gerade - straight - recto

18,5 cm / 7 1/4"

HH 500-02 2 mm

HH 500-04 4 mm

HH 500-07 7 mm

HH 500-09 9 mm

HH 500-12 12 mm



### Fomon

Osteotome - Osteotomes - Osteótomos  
gerade - straight - recto

16 cm / 6 1/4"

HH 620-04 4 mm

HH 620-05 5 mm

HH 620-06 6 mm

HH 620-07 7 mm



### Cinelli

Osteotome - Osteotomes - Osteótomos  
gerade - straight - recto

16 cm / 6 1/4"

HH 630-10 10 mm

HH 630-12 12 mm


HH 630-14 14 mm

HH 630-16 16 mm



**Nasenraspeln**  
**Nasal Rasps**  
**Raspas nasales**




**Maltz**  
**HH 720-18** 18 cm / 7"



**Fomon**  
 doppelendig - double-ended - dos extremos  
 fein/grob - fine/coarse - fina/gruesa  
**HH 730-21** 21 cm / 8 1/4"




**Aufricht**  
 auf Schub schneidend  
 pushing cut  
 corte empujando  
**HH 751-21** 21 cm / 8 1/4"



**HH 751-21**

auf Zug schneidend  
 pulling cut  
 corte tirando  
**HH 753-21** 21 cm / 8 1/4"



**HH 753-21**

**Nasenraspeln**  
**Nasal Rasps**  
**Raspas nasales**



**GEOMED GOLD-CUT<sup>®</sup>**

**Fomon**

Hartmetallraspel - TC Rasp - TC raspa  
doppelseitig - double-ended - dos extremos  
20 cm / 8"


**HH 770-01** # 1 / 2

**HH 770-02** # 3 / 4

**HH 770-03** # 5 / 6


**HH 770-04** # 7 / 8

**HH 770-05** # 9 / 10




# 1

# 2




# 3

# 4




# 5

# 6



# 7

# 8



# 9

# 10



**Lewis**

auf Zug schneidend  
pulling cut  
corte tirando

fein - fine - fina

**HH 781-18** 18 cm / 7"



**HH 781-18**



**HH 783-18**

grob - coarse - gruesa

**HH 783-18** 18 cm / 7"

## Nasensägen, Antrumspülkatheter

### Nasal Saws, Antrum Irrigation Cannulas

### Sierras nasales, canulas para irrigación de los senos maxilares



**Joseph**

gerade - straight - recta

**HH 850-01** 18 cm / 7"



**HH 850-01**

bajonettförmig - bayonet shaped - forma bayoneta

rechts gebogen - curved to right - curva a derecha

**HH 851-01** 18 cm / 7"



**HH 851-01**

bajonettförmig - bayonet shaped - forma bayoneta

links gebogen - curved to left - curva a izquierda

**HH 851-02** 18 cm / 7"



**HH 851-02**



**Kilian**

kurze Biegung - short curved - curva corta

15 cm / 6"

**HH 981-15** ø 1,5 mm

**HH 981-25** ø 2,5 mm

**HH 981-30** ø 3,0 mm

**HH 981-40** ø 4,0 mm



**Kilian**

lange Biegung - long curved - curva larga

14,5 cm / 5¾"


**HH 983-15** ø 1,5 mm

**HH 983-25** ø 2,5 mm

**HH 983-30** ø 3,0 mm

**HH 983-40** ø 4,0 mm

**Nasenpolypenzangen, Siebbeinzangen**  
**Nasal Polypus Forceps, Ethmoidal Forceps**  
**Pinzas para pólipos nasales, pinzas para etmoides**



**Hartmann**  
**HK 101-16** 16 cm / 6¼"  
**HK 101-20** 20 cm / 8"



**Weil-Blakesley**  
ohne Hals - without neck - sin cuello

Schaft - shaft - tubo 12,5 cm / 4¾"

**HK 370-01**  
**HK 370-02**  
**HK 370-03**  
**HK 370-04**  
**HK 370-05**

aufwärts - up - arriba 45°

**HK 371-01**  
**HK 371-02**  
**HK 371-03**  
**HK 371-04**  
**HK 371-05**



**HK 370-01**  
**HK 371-01**  
**HK 376-01**  
**HK 377-01**



**HK 370-02**  
**HK 371-02**  
**HK 376-02**  
**HK 377-02**



**HK 370-03**  
**HK 371-03**



**HK 370-04**  
**HK 371-04**



**HK 370-05**  
**HK 371-05**  
**HK 376-05**  
**HK 377-05**

**Weil-Blakesley**

durchschneidend - through-cut - cortante  
gerade - straight - recta

Schaft - shaft - tubo 10,5 cm / 4¼"

**HK 376-01**  
**HK 376-02**  
**HK 376-05**


aufwärts - up - arriba 45°

**HK 377-01**  
**HK 377-02**  
**HK 377-05**

**Nasenconchotome, Kieferhöhlenstanzen**

**Nasal Cutting Forceps, Antrum Punches**

**Pinzas cortantes nasales, pinzas cortantes para senos maxilares**



**Hartmann**

19 cm / 7½"

Schaft - shaft - tubo 12 cm / 4¾"

**HK 400-05** ø 5 mm

**HK 400-07** ø 7 mm

**HK 400-09** ø 9 mm

**HK 400-11** ø 11 mm



**HK 400-05**




**HK 400-07**



**HK 400-09**



**HK 400-11**



**Hajek-Kofler**

Schaft - shaft - tubo 14 cm / 5½"

aufwärts - up - arriba


**HK 470-03** 3,5 mm x 3,5 mm

**HK 470-04** 4 mm x 4 mm

abwärts - down - abajo

**HK 472-03** 3,5 mm x 3,5 mm


**HK 472-04** 4 mm x 4 mm




**HK 470-03**



**HK 470-04**



**HK 472-03**



**HK 472-04**



**Septumrichtzangen**  
**Septum Straightening Forceps**  
**Pinzas para enderezar el septo**



**Walsham-Universal**

flache Blätter

flat blades

hojas planas

**HK 800-23** 23 cm / 9"



**Walsham**

linkes Blatt geformt

left blade grooved

hoja izquierda formada

**HK 801-23** 23 cm / 9"

rechtes Blatt geformt

right blade grooved

hoja derecha formada

**HK 803-23** 23 cm / 9"



**Ash**



aufgewinkelt

angled on flat

angulada de plano

**HK 811-23** 23 cm / 9"


**Septumzangen und -scheren**  
**Septum Forceps and Scissors**  
**Pinzas y tijeras para el septo**



**Jansen-Middleton**


Septumzange, durchschneidend  
Septum forceps, through-cutting  
Pinza para el septo, corte perforando  
**HK 821-19** 19 cm / 7½"

Septumzange, löffelförmiges Maul  
Septum forceps, spoon-shaped jaw  
Pinza para el septo, bocas forma cuchara  
**HK 823-19** 19 cm / 7½"



**Caplan**

Septumschere, gezahnt  
Septum scissors, serrated  
Tijera para el septo, dentada  
**HK 831-20** 20 cm / 8"



**Rubin**

Knorpelquetschzange, mit Schutzkappe  
Cartilage morselising forceps, with cap  
Pinza aplastanda para cartílagos, con protección  
**HK 851-21** 21 cm / 8¼"



Schutzkappe - protection cap - protección  
**HK 852-01**

**Ohrtrichter**  
**Ear Specula**  
**Espéculos auriculares**



**HO 100-03**



**HO 100-04**



**HO 100-05**



**HO 100-06**



**HO 100-07**

**Hartmann**

rund - round - redondo

**HO 100-03** # 0 ø 3,0 mm

**HO 100-04** # 1 ø 4,5 mm

**HO 100-05** # 2 ø 5,5 mm

**HO 100-06** # 3 ø 6,5 mm

**HO 100-07** # 4 ø 7,5 mm

Satz - set - juego

für Kinder - for children - para niños

**HO 101-02**

bestehend aus - consisting of - compuesto de:

**HO 100-03**

**HO 100-04**

**HO 100-05**

Satz - set - juego

für Erwachsene - for adults - para adultos

**HO 101-03**

bestehend aus - consisting of - compuesto de:

**HO 100-05**

**HO 100-06**

**HO 100-07**

Satz - set - juego

für Erwachsene - for adults - para adultos

**HO 101-04**

bestehend aus - consisting of - compuesto de:

**HO 100-04**

**HO 100-05**

**HO 100-06**

**HO 100-07**

**Stimmgabeln**  
**Tuning Forks**  
**Diapasones**



**Hartmann**  
**HO 218-10**    c    128



**HO 218-11**    c<sup>1</sup>    256



**Hartmann**  
**HO 220-12**    c<sup>2</sup>    512



**HO 220-13**    c<sup>3</sup>    1024



**HO 220-14**    c<sup>4</sup>    2048



**HO 220-15**    c<sup>5</sup>    4096



**Hartmann**  
**HO 221-00**  
 Satz komplett in Etui  
 Set complete in case  
 Juego completo en estuche

bestehend aus - including - compuesto de:

- HO 218-10**    c    128
- HO 218-11**    c<sup>1</sup>    256
- HO 220-12**    c<sup>2</sup>    512
- HO 220-13**    c<sup>3</sup>    1024
- HO 220-14**    c<sup>4</sup>    2048



**Rydel-Seiffer**  
 Dämpfer abnehmbar, mit Fuß  
 Detachable clamps, with foot  
 Sordinas demontables, con pie  
**HO 230-01**    c 64 / c 128

**Ohrkatheter und -spritzen**  
**Eustachian Catheters, Ear Syringes**  
**Tubos y jeringas auriculares**



**Hartmann**

14,5 cm / 5¾"

**HO 301-20**    ø 2,0 mm

**HO 301-25**    ø 2,5 mm

**HO 301-30**    ø 3,0 mm

**HO 301-35**    ø 3,5 mm

**HO 301-40**    ø 4,0 mm

**HO 301-50**    ø 5,0 mm



**Itard**

14,5 cm / 5¾"


**HO 311-30**    ø 3,0 mm

**HO 311-35**    ø 3,5 mm

**HO 311-40**    ø 4,0 mm

**HO 311-45**    ø 4,5 mm

**HO 311-50**    ø 5,0 mm



Ohrspritze,  
Ganzmetall, vernickelt, Luer-Lock-Ansatz,  
mit 2 Kanülen und 1 Schutzteller.

Ear syringe,  
all metal, chromium-plated, Luer-Lock cone,  
with 2 cannulas and 1 protection shield.

Jeringa auricular,  
de metal, cromada, conexión Luer-Lock,  
con 2 cánulas y 1 protector.




komplett - complete - completo

**HO 350-05**    50 ml

**HO 350-10**    100 ml

**Ohrpolypenschnürer**  
**Ear Polypus Snares**  
**Cierra-nudos auriculares**



Schnürdraht, Ring mit 10 m  
 Snare wire, coil of 10 m  
 Alambre, rollo de 10 m  
**HO 405-03** ø 0,3 mm  
**HO 405-04** ø 0,4 mm




Schnürdraht, geformt, Dtzd.  
 Snare wire, shaped, doz.  
 Alambre, formado, dc.  
**HO 407-03** ø 0,3 mm



**Krause**  
**HO 411-16** 16 cm / 6¼"



**Krause-Voss**  
**HO 421-16** 16 cm / 6¼"



**Lange-Wilde**  
**HO 431-16** 16 cm / 6¼"


## Ohrküretten und -schlingen Ear Curettes and Loops Curetas y asas auriculares



### Buck

16,5 cm / 6½"

gerade - straight - recta



stumpf - blunt - roma

HO 700-01	# 000
HO 700-02	# 00
HO 700-03	# 0
HO 700-04	# 1
HO 700-05	# 2
HO 700-06	# 3
HO 700-07	# 4




scharf - sharp - cortante

HO 706-01	# 000
HO 706-02	# 00
HO 706-03	# 0
HO 706-04	# 1
HO 706-05	# 2
HO 706-06	# 3
HO 706-07	# 4



### Billeau

16 cm / 6¼"




HO 710-01	# 1
HO 710-02	# 2
HO 710-03	# 3



### Shapleigh

16 cm / 6¼"



gerade - straight - recta

HO 740-01	# 1
HO 740-02	# 2

gebogen - curved - curva

HO 741-01	# 1
HO 741-02	# 2

## Parazentesenadeln

## Paracentesis Needles

## Agujas para paracentesis



### Politzer


gerade - straight - recta

vertikaler Schnitt

vertical cut

corte vertical

**HO 770-16** 16 cm / 6¼"



### Politzer


abgewinkelt - angled - angulada

vertikaler Schnitt

vertical cut

corte vertical

**HO 771-16** 16 cm / 6¼"



horizontaler Schnitt

horizontal cut

corte horizontal

**HO 773-16** 16 cm / 6¼"



### Lucae

bajonettförmig - bayonet shaped - forma bayoneta

vertikaler Schnitt

vertical cut

corte vertical

**HO 783-18** 18 cm / 7"



horizontaler Schnitt

horizontal cut

corte horizontal

**HO 785-18** 18 cm / 7"



**Ohrhäkchen**  
**Ear Hooks**  
**Ganchos auriculares**



**Wagener**

14 cm / 5½"

- HO 815-01** # 1
- HO 815-02** # 2
- HO 815-03** # 3
- HO 815-04** # 4
- HO 815-05** # 5



**HO 815-01**



**HO 815-02**



**HO 815-03**



**HO 815-04**



**HO 815-05**



**Day**

16 cm / 6¼"

- HO 825-01** # 1
- HO 825-02** # 2



**HO 825-01**



**HO 825-02**




**Lucae**

14 cm / 5½"

- HO 855-01** # 1
- HO 855-02** # 2
- HO 855-03** # 3



**HO 855-01**



**HO 855-02**



**HO 855-03**




**Groß**

Fremdkörperinstrument

Foreign body spud

Instrumento para cuerpos extraños



stumpf - blunt - romo

**HO 861-12** 12 cm / 4¾"

**Ohrpolypenzängchen, Fremdkörperzängchen**

**Ear Polypus Forceps, Foreign Body Forceps**

**Pinzas para pólipos auriculares, pinzas para cuerpos extraños**



HS 101-13

HS 111-13

**Hartmann**

feines Modell - fine pattern - modelo delgado


HS 101-13 13 cm / 5"

Standardmodell

Standard pattern


Modelo standard

HS 111-13 13 cm / 5"



**Buck**

HS 120-13 13,5 cm / 5¼"



**Quire**

HS 130-11 11,5 cm / 4½"

**Ohrzängchen**  
**Ear Forceps**  
**Pinzas auriculares**



Längenangabe = Schaftlänge

Length indication = length of shaft

Indicación de longitud = longitud de tubo




**Hartmann-Wullstein**

**HS 200-08** 8,5 cm / 3¼"



ø 1 mm

**HS 210-08** 8,5 cm / 3¼"



**Hartmann**

**HS 220-08** 8,5 cm / 3¼"

**HS 220-14** 14 cm / 5½"



**HS 230-08** 8,5 cm / 3¼"

**HS 230-14** 14 cm / 5½"



**HS 240-08**

8,5 cm / 3¼"



ø 2,0 mm

**HS 250-20** 8,5 cm / 3¼"

ø 2,5 mm

**HS 250-25** 8,5 cm / 3¼"

ø 3,0 mm

**HS 250-30** 8,5 cm / 3¼"



**Struempel**

**HS 270-08** 8,5 cm / 3¼"



**Hartmann-Hoffmann**

ø 2,0 mm

**HS 280-08** 8,5 cm / 3¼"

## Knochenstanzen

### Bone Punches

### Punzonadoras de huesos

*Längenangabe = Schaftlänge*

*Length indication = length of shaft*

*Indicación de longitud = longitud de tubo*



#### **Citelli**

7,5 cm / 3"

Breite - width - anchura 1,5 mm  
**HS 620-15**

Breite - width - anchura 2,2 mm  
**HS 620-22**

Breite - width - anchura 3 mm  
**HS 620-30**



#### **Kerrison**

9 cm / 3½"

Breite - width - anchura 3 mm  
**HS 640-03**

Breite - width - anchura 4 mm  
**HS 640-04**

Breite - width - anchura 5 mm  
**HS 640-05**

Breite - width - anchura 6 mm  
**HS 640-06**

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

## ZH 101-00

### HNO Instrumenten-Set

### ENT Instrument Set

### Set Instrumental Otorrinolaringología

AA 131-12	2	WAGNER SURG SCISSOR BL/BL CVD 12CM
AB 310-11	2	STEVENS DISSECT SCISSOR BL/BL STR 11CM
AB 311-11	2	STEVENS DISSECT SCISSOR BL/BL CVD 11CM
AB 341-17	2	PRINCE TONSIL SCISS CVD 17.5CM
AB 957-24	1	BRUNS PLASTER SHEAR SERRATED 24CM
AE 486-12	2	ADSON TISSUE FORCEPS 1X2T STR 12CM
AF 515-20	2	GRUENWALD NASAL DRESSING FORCEPS 20CM
AF 535-14	1	LUCAE EAR FORCEPS 14CM
AF 551-18	2	TROELTSCH (WILDE) EAR FORCEPS 18CM
AI 100-12	12	HALSTEAD-MOSQUITO HEMOST FCPS STR 12.5CM
AI 101-12	12	HALSTEAD-MOSQUITO HEMOST FCPS CVD 12.5CM
AI 130-14	8	CRILE HEMOSTATIC FORCEPS STR 14CM
AI 131-14	8	CRILE HEMOSTATIC FORCEPS CVD 14CM
AI 333-16	8	PEAN (ROCHESTER) HEMOSTAT FCPS CVD 16CM
AK 151-11	24	BACKHAUS TOWEL FORCEPS 11CM
AL 200-25	6	FOERSTER SPONGE HOLD FCPS SERR STR 25CM
AR 220-15	1	GOLD-GRIP CRILE-WOOD NDHL P04 15CM
AT 121-08	2	FRAZIER SUCTION TUBE 2.7MM 19CM
AT 121-12	1	FRAZIER SUCTION TUBE 4.0MM 19CM
AT 200-00	1	YANKAUER SUCTION TUBE COMPLETE 27CM
AU 311-16	1	COAKLEY TROCAR DIA 2.25MM 16CM
AW 110-14	2	BUTTON-END PROBE DIA 2.0MM 14.5CM
AW 150-14	2	GROOVED DIRECTOR 14.5CM
AW 300-16	2	FARRELL COTTON APPLIC RD SERR 1.1MM/16CM
AW 300-14	2	FARRELL COTTON APPLIC RD SERR 0.9MM/14CM
AX 044-16	2	GRAEFE HOOK SHARP 16CM
AX 240-14	2	CUSHING RETRACTOR 14MM 20CM
AX 630-00	1	PARKER-LANGENBECK RETRACT SET #1-2 21CM
AZ 120-16	1	WEITLANER RETRACTOR 3X4T BLUNT 16.5CM
MC 110-15	8	ALLIS INTESTINAL FORCEPS 4X5T 15CM
MC 350-16	4	BABCOCK ATRAUMAT INTESTINAL FORCEPS 16CM
HA 105-12	2	LARYNGEAL MIRROR W/O HANDLE #1 DIAM 12MM
HA 105-16	2	LARYNGEAL MIRROR W/O HANDLE #3 DIAM 16MM
HG 160-14	2	WIEN NASAL SPECULA 30MM 14.5CM
HG 200-01	2	HARTMANN-HALLE NASAL SPECULA 25MM 15CM
HH 201-01	1	COAKLEY ANTRUM CURETTE #1 17CM
HH 201-02	1	COAKLEY ANTRUM CURETTE #2 17CM
HH 201-03	1	COAKLEY ANTRUM CURETTE #3 17CM
HH 201-04	1	COAKLEY ANTRUM CURETTE #4 17CM
HH 201-05	1	COAKLEY ANTRUM CURETTE #5 17CM
HH 201-06	1	COAKLEY ANTRUM CURETTE #6 17CM
HH 315-18	1	FREER SEPTUM ELEVAT 5MM/5MM 18CM
HH 390-22	1	KILLIAN SEPTUM ELEVATOR
HH 751-21	1	AUFRIECHT NASAL RASP PUSHING CUT 21CM
HH 753-21	1	AUFRIECHT NASAL RASP PULLING CUT 21CM
HH 971-30	2	EICKEN ANTRUM IRRIG CANNULA 3.0MM 14.5CM
HK 821-19	1	JANSEN-MIDDLETON SEPTUM FORCEPS 19CM
HK 823-19	1	JANSEN-MIDDLETON SEPTUM FORCEPS 19CM
HO 218-10	1	HARTMANN TUNING FORK C 128HZ
HO 218-11	1	HARTMANN TUNING FORK C1 256HZ
HO 411-16	1	KRAUSE EAR POLYPUS SNARE 16CM

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

HO 405-03	10	SNARE WIRE COIL OF 10M DIA 0.3MM
HO 861-12	2	GROSS FOREIGN BODY SPUD BLUNT 12CM
HS 111-14	2	HARTMANN EAR POLYPUS FCPS STAND PAT 14CM
HS 120-13	1	BUCK FOREIGN BODY FORCEPS 13.5CM
HS 200-08	2	HARTMANN-WULLSTEIN EAR FORCEPS 8,5CM
HS 640-03	1	KERRISON BONE PUNCH 3MM 9CM
HS 640-04	1	KERRISON BONE PUNCH 4MM 9CM
HS 640-05	1	KERRISON BONE PUNCH 5MM 9CM
OD 142-03	1	PARTSCH CHISEL 3MM 17CM
OD 142-04	1	PARTSCH CHISEL 4MM 17CM
OD 142-06	1	PARTSCH CHISEL 6MM 17CM
OD 144-03	1	PARTSCH GOUGE 3MM 17CM
OD 144-04	1	PARTSCH GOUGE 4MM 17CM
OD 144-06	1	PARTSCH GOUGE 6MM 17CM
OD 580-02	1	SPRATT BONE CURETTE OVAL STR #000 17CM
OD 580-04	1	SPRATT BONE CURETTE OVAL STR #0 17CM
OD 580-06	1	SPRATT BONE CURETTE OVAL STR #2 17CM
OD 580-08	1	SPRATT BONE CURETTE OVAL STR #4 17CM
OR 121-14	1	FRIEDMAN BONE RONGEUR 30DEG 14CM

## ZH 102-00

### HNO-Untersuchungs-Set

### ENT Examination Set

### Set de Examen Otorrinolaringología

AF 551-18	1	TROELTSCH (WILDE) EAR FORCEPS 18CM
AW 110-14	1	BUTTON-END PROBE DIA 2.0MM 14.5CM
AW 300-14	1	FARRELL COTTON APPLIC RD SERR 0.9MM/14CM
AW 300-18	1	FARRELL COTTON APPLIC RD SERR 1.2MM/18CM
HA 105-16	1	LARYNGEAL MIRROR W/O HANDLE #3 DIAM 16MM
HA 105-18	1	LARYNGEAL MIRROR W/O HANDLE #4 DIAM 18MM
HA 105-20	1	LARYNGEAL MIRROR W/O HANDLE #5 DIAM 20MM
HA 110-10	1	HANDLE FOR LARYNGEAL MIRRORS 10CM
HC 070-19	1	BRUENINGS TONGUE DEPRESSOR 19CM
HC 101-14	1	TOBOLD TONGUE DEPRESSOR 14.5CM
HG 160-14	1	WIEN NASAL SPECULA 30MM 14.5CM
HK 101-20	1	HARTMANN NASAL POLYPUS FORCEPS 20CM
HO 131-03	1	TOYNBEE SET OF 3 EAR SPECULA #2-4
HO 770-16	1	POLITZER PARACENT NEEDLE STR VERT 16CM
HO 861-12	1	GROSS FOREIGN BODY SPUD BLUNT 12CM
ID 660-00	1	COMBINATION DIAGONSTIC SET COMPLET

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

### ZH 103-00

#### Instrumenten-Set für die Trommelfellplastik Myringoplasty Instrument Set Set Extracción de Cuerpo Extraño

AF 515-20	1	GRUENWALD NASAL DRESSING FORCEPS 20CM
AS 100-03	1	SCALPEL HANDLE NO 3 STANDARD
AS 100-07	1	SCALPEL HANDLE NO 7
AT 121-06	1	FRAZIER SUCTION TUBE 2.0MM 19CM
AT 121-08	1	FRAZIER SUCTION TUBE 2.7MM 19CM
AW 300-16	1	FARRELL COTTON APPLIC RD SERR 1.1MM/16CM
HO 100-05	1	HARTMANN EAR SPECULA ROUND #2 5.5MM
HO 100-06	1	HARTMANN EAR SPECULA ROUND #3 6.5MM
HO 710-02	1	BILLEAU EAR CURETTE #2 16CM
HO 861-12	1	GROSS FOREIGN BODY SPUD BLUNT 12CM
HS 200-08	1	HARTMANN-WULLSTEIN EAR FORCEPS 8,5CM

### ZH 104-00

#### Instrumenten-Set für die Polypenexzision des Gehörgangs Aural Polypectomy Set Set Polipectomía

AF 515-20	1	GRUENWALD NASAL DRESSING FORCEPS 20CM
AS 100-03	1	SCALPEL HANDLE NO 3 STANDARD
AS 100-07	1	SCALPEL HANDLE NO 7
AT 121-06	1	FRAZIER SUCTION TUBE 2.0MM 19CM
AT 121-08	1	FRAZIER SUCTION TUBE 2.7MM 19CM
AW 300-16	1	FARRELL COTTON APPLIC RD SERR 1.1MM/16CM
HO 100-05	1	HARTMANN EAR SPECULA ROUND #2 5.5MM
HO 100-06	1	HARTMANN EAR SPECULA ROUND #3 6.5MM
HO 100-07	1	HARTMANN EAR SPECULA ROUND #4 7.5MM
HO 411-16	1	KRAUSE EAR POLYPUS SNARE 16CM
HO 710-02	1	BILLEAU EAR CURETTE #2 16CM
HO 861-12	1	GROSS FOREIGN BODY SPUD BLUNT 12CM
HS 200-08	1	HARTMANN-WULLSTEIN EAR FORCEPS 8,5CM

### ZH 105-00

#### Instrumenten-Set für die Fremdkörperentfernung im Ohr Ear Foreign Body Remove Set Set Remoción Cuerpo Extraño del Oído

AF 535-14	1	LUCAE EAR FORCEPS 14CM
HO 101-04	1	HARTMANN EAR SPECULA SET FOR ADULTS
HO 710-01	1	BILLEAU EAR CURETTE #1 16CM
HO 710-02	1	BILLEAU EAR CURETTE #2 16CM
HS 200-08	1	HARTMANN-WULLSTEIN EAR FORCEPS 8,5CM

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

### ZH 107-00

#### Mastoidektomie-Set Mastoidectomy Set Set de Mastoidectomía

AA 170-14	1	MAYO DISSECTING SCISSOR STR 14.5CM
AE 486-12	2	ADSON TISSUE FORCEPS 1X2T STR 12CM
AI 332-16	2	PEAN (ROCHESTER) HEMOSTAT FCPS STR 16CM
AI 100-12	3	HALSTEAD-MOSQUITO HEMOST FCPS STR 12.5CM
AL 200-18	1	FOERSTER SPONGE HOLD FCPS SERR STR 18CM
AR 120-13	1	GOLD-GRIP HALSEY NEEDLE HOLDER P04 13CM
AS 100-03	1	SCALPEL HANDLE NO 3 STANDARD
AS 100-07	1	SCALPEL HANDLE NO 7
AT 121-08	1	FRAZIER SUCTION TUBE 2.7MM 19CM
AT 121-10	1	FRAZIER SUCTION TUBE 3.3MM 19CM
AW 110-14	1	BUTTON-END PROBE DIA 2.0MM 14.5CM
AW 300-16	2	FARRELL COTTON APPLIC RD SERR 1.1MM/16CM
AX 044-16	1	GRAEFE HOOK SHARP 16CM
AX 320-03	2	VOLKMANN RETRACTOR 3 PRONGS BLUNT 21.5CM
AX 660-00	1	FARABEUF RETRACTOR SET #1-2 15CM
AX 680-03	2	RICHARDSON RETRACTOR 38X38MM 24CM
DI 200-14	1	STERNBERG LIP RETRACTOR
HC 080-14	1	WIEDER TONGUE DEPRESSOR 14CM
HC 131-12	1	DENHART MOUTH GAG 12CM
HG 160-14	1	WIEN NASAL SPECULA 30MM 14.5CM
HG 651-20	1	COTTLE-NEIVERT HOOK
HH 201-01	1	COAKLEY ANTRUM CURETTE #1 17CM
HH 201-02	1	COAKLEY ANTRUM CURETTE #2 17CM
HH 201-03	1	COAKLEY ANTRUM CURETTE #3 17CM
HH 250-23	1	COTTLE SEPTUM ELEVAT 4.5MM/4.5MM 23CM
HH 241-16	1	JOSEPH PERIOSTEAL ELEVATOR 16CM
HH 333-16	1	FREER SEPTUM CHISEL
HH 500-04	1	COTTLE SEPTUM OSTEOTOME STR 4MM 18.5CM
HH 500-07	1	COTTLE SEPTUM OSTEOTOME STR 7MM 18.5CM
HS 220-16	1	HARTMANN EAR POLYPUS FCPS 16CM
NS 520-15	1	KERRISON LAMIN PUNCH UP 3MM 15CM
MC 120-19	2	ALLIS INTESTINAL FORCEPS 5X6T 19CM
OD 920-17	1	MEAD Mallet W PLASTIC FACES 320G 17CM
OR 431-18	1	BEYER BONE RONGEUR CVD 18CM

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

### ZH 108-00

#### Tonsillektomie- und Anedoidektomie-Set Tonsillectomy / Anoidectomy Set Set de Adenoamigdalectomía y Adenoidectomía

AD 551-18	2	GOLD-CUT METZENB SCISSOR BL/BL CVD 18CM
AI 425-18	1	SCHNIDT DISSECTING CLAMP STRONG CVD 18CM
AK 151-08	4	BACKHAUS TOWEL FORCEPS 8CM
AL 200-18	2	FOERSTER SPONGE HOLD FCPS SERR STR 18CM
AL 201-18	2	FOERSTER SPONGE HOLD FCPS SERR CVD 18CM
AR 120-13	2	GOLD-GRIP HALSEY NEEDLE HOLDER P04 13CM
AS 100-07	1	SCALPEL HANDLE NO 7
AT 200-00	1	YANKAUER SUCTION TUBE COMPLETE 27CM
HC 170-09	1	JENNINGS MOUTH GAG 9CM
HD 171-22	1	HURD DISSECTOR 22.5CM
HD 431-18	1	WHITE TONSIL FORCEPS CVD 18CM
HD 495-30	1	SNARE WIRE SHAPED DOZEN 30MMX100MM
HD 500-01	1	TYDING TONSIL SNARE COMPL W STRAIGHT TIP
HD 530-13	1	HARTMANN TONSIL PUNCH DIAM 13MM 14CM
HD 561-02	1	BARNHILL ADENOID CURETTE STR #1 22CM
HD 561-03	1	BARNHILL ADENOID CURETTE STR #2 22CM
HD 561-04	1	BARNHILL ADENOID CURETTE STR #3 22CM
HD 650-00	1	SLUDER-BALLENGER TONSILLECTOME COMP
NC 110-19	1	LOVE NERVE ROOT RETRACTOR STR 19CM
NC 111-19	1	LOVE NERVE ROOT RETRACTOR ANGLED 19CM

### ZH 109-00

#### Instrumenten-Set für die Gaumenplastik Instrument Set for Palatoplasty Set de Paladar

AI 425-18	2	SCHNIDT DISSECTING CLAMP STRONG CVD 18CM
AR 200-20	1	GOLD-GRIP MAYO-HEGAR NDHL P05 20CM
AS 100-07	1	SCALPEL HANDLE NO 7
HC 101-14	1	TOBOLD TONGUE DEPRESSOR 14.5CM
HC 131-12	1	DENHART MOUTH GAG 12CM
HD 171-22	2	HURD DISSECTOR 22.5CM
HH 500-07	1	COTTLE SEPTUM OSTEOTOME STR 7MM 18.5CM
NA 030-00	1	BUNNELL HAND-DRILL SET COMPLETE
NC 110-19	1	LOVE NERVE ROOT RETRACTOR STR 19CM
OD 451-06	1	WILLIGER ELEVATOR LIGHT CURVE 6.5MM 17CM

### ZH 110-00

#### Instrumenten-Set für die Kehlkopfspaltung Laryngofissure Set Set Laringofisura

AE 160-12	2	ADSON DRESSING FORCEPS STR 12CM
AE 250-21	1	POTTS-SMITH DRESSING FCPS STR 21CM
AE 486-12	2	ADSON TISSUE FORCEPS 1X2T STR 12CM
AT 121-08	1	FRAZIER SUCTION TUBE 2.7MM 19CM
HD 171-22	1	HURD DISSECTOR 22.5CM
HH 315-18	1	FREER SEPTUM ELEVAT 5MM/5MM 18CM

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

### ZH 111-00

#### Tracheostomie-Set Tracheostomy Set Set de Traqueostomía

AA 102-14	1	SURGICAL SCISSOR SH/BL STR 14.5CM
AD 511-17	1	GOLD-CUT MAYO DISSECT SCISS CVD 17CM
AD 551-18	1	GOLD-CUT METZENB SCISSOR BL/BL CVD 18CM
AD 601-14	1	GOLD-CUT METZENBAUM SCISSOR CVD 14.5CM
AE 160-12	1	ADSON DRESSING FORCEPS STR 12CM
AE 406-14	1	TISSUE FORCEPS 1X2T STR 14.5CM
AE 406-14	1	TISSUE FORCEPS 1X2T STR 14.5CM
AE 486-12	1	ADSON TISSUE FORCEPS 1X2T STR 12CM
AI 100-12	2	HALSTEAD-MOSQUITO HEMOST FCPS STR 12.5CM
AI 101-12	3	HALSTEAD-MOSQUITO HEMOST FCPS CVD 12.5CM
AI 131-14	4	CRILE HEMOSTATIC FORCEPS CVD 14CM
AK 151-13	4	BACKHAUS TOWEL FORCEPS 13CM
AL 200-25	4	FOERSTER SPONGE HOLD FCPS SERR STR 25CM
AR 200-18	1	GOLD-GRIP MAYO-HEGAR NDHL P05 18CM
AS 100-03	1	SCALPEL HANDLE NO 3 STANDARD
AT 121-08	1	FRAZIER SUCTION TUBE 2.7MM 19CM
AT 200-00	1	YANKAUER SUCTION TUBE COMPLETE 27CM
AX 100-02	2	HOOK FINE PATTERN 2 PRONGS BLUNT 16CM
AX 110-02	2	STANDARD HOOK 2 PRONGS BLUNT 16CM
AX 154-03	2	SENN-MILLER RETRACTOR SHARP 20X6MM 16CM
HA 121-14	1	LABORDE TRACHEAL DILATOR 14CM
HA 214-15	2	JACKSON TRACHEAL HOOK 15CM
HA 425-05	1	CHEV-JACKSON TRACH TUBE W OBTUR #1 5MM
HA 425-06	1	CHEV-JACKSON TRACH TUBE W OBTUR #2 6MM
HA 425-10	1	CHEV-JACKSON TRACH TUBE W OBTUR #6 10MM
HA 425-11	1	CHEV-JACKSON TRACH TUBE W OBTUR #7 11MM
MC 120-15	2	ALLIS INTESTINAL FORCEPS 5X6T 15CM

### ZH 112-00

#### Thyroidektomie-Set Thyroidectomy Set Set de Tiroidectomía

AI 130-14	24	CRILE HEMOSTATIC FORCEPS STR 14CM
AI 131-14	24	CRILE HEMOSTATIC FORCEPS CVD 14CM
AI 611-19	6	LAHEY GALL DUCT CLAMP CVD 19CM
AN 510-21	1	DESCHAMPS LIGAT NEEDLE RIGHT CVD BL 21CM
AN 511-21	1	DESCHAMPS LIGAT NEEDLE LEFT CVD BL 21CM
AX 240-18	2	CUSHING RETRACTOR 18MM 20CM
AX 270-20	2	LAHEY RETRACTOR 20 X 7 MM 20CM

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

Artikel-Nr.	Menge quantity cantidad	Beschreibung Description descripción
-------------	-------------------------------	--

## ZH 114-00

### Instrumenten-Set für die Nasenplastik

#### Rhinoplasty Set

#### Set de Rinoseptoplastía

AA 102-14	1	SURGICAL SCISSOR SH/BL STR 14.5CM
AA 201-18	1	METZENBAUM DISS SCISSOR BL/BL CVD 18CM
AB 111-18	1	HEYMANN NASAL SCISSOR 18CM
AB 310-11	1	STEVENS DISSECT SCISSOR BL/BL STR 11CM
AE 100-13	1	DRESSING FORCEPS STANDARD STR 13CM
AE 160-12	1	ADSON DRESSING FORCEPS STR 12CM
AE 486-12	1	ADSON TISSUE FORCEPS 1X2T STR 12CM
AE 536-12	1	ADSON-BROWN TISSUE FORCEPS 9X9T 12CM
AF 515-20	1	GRUENWALD NASAL DRESSING FORCEPS 20CM
AI 100-12	3	HALSTEAD-MOSQUITO HEMOST FCPS STR
AI 101-12	3	HALSTEAD-MOSQUITO HEMOST FCPS CVD
AI 120-15	1	LERICHE HEMOSTATIC FORCEPS STR 15CM
AK 151-08	4	BACKHAUS TOWEL FORCEPS 8CM
AK 151-11	4	BACKHAUS TOWEL FORCEPS 11CM
AL 200-18	2	FOERSTER SPONGE HOLD FCPS SERR STR 18CM
AN 300-14	1	MAYO-HEGAR NEEDLE HOLDER 14CM
AN 300-16	1	MAYO-HEGAR NEEDLE HOLDER 16CM
AN 300-18	1	MAYO-HEGAR NEEDLE HOLDER 18CM
AS 100-03	1	SCALPEL HANDLE NO 3 STANDARD
AS 100-04	1	SCALPEL HANDLE NO 4
AT 121-07	1	FRAZIER SUCTION TUBE 7 FR
AT 121-09	1	FRAZIER SUCTION TUBE 9 FR
AT 200-00	1	YANKAUER SUCTION TUBE COMPLETE 27CM
AX 154-03	2	SENN-MILLER RETRACTOR SHARP 20X6MM 16CM
AX 665-00	1	FARABEUF RETRACTOR SET #1-2 12CM
HG 200-01	1	HARTMANN-HALLE NASAL SPECULA 25MM 15CM
HG 200-02	1	HARTMANN-HALLE NASAL SPECULA 30MM 15CM
HG 200-03	1	HARTMANN-HALLE NASAL SPECULA 35MM 15CM
HH 301-19	1	FREER SEPTUM KNIFE
HH 315-18	1	FREER SEPTUM ELEVAT 5MM/5MM 18CM
HH 315-18	2	FREER SEPTUM ELEVAT 5MM/5MM 18CM
HH 390-22	1	KILLIAN SEPTUM ELEVATOR
HH 501-03	1	COTTLE SEPTUM OSTEOTOME CVD 18X3MM
HH 751-21	1	AUFRIEHT NASAL RASP PUSHING CUT 21CM
HH 851-02	1	JOSEPH NASAL SAW BAYONET LEFT CV 18CM
HK 260-02	1	LUC SEPTUM FCPS FIG 2
HK 801-23	1	WALSHAM SEPT FCPS LEFT BLADE GROOV 23CM
HK 803-23	1	WALSHAM SEPT FCPS RIGHT BLADE GROOV
HK 811-23	1	ASH SEPTUM STRAIGHTENING FCPS ANGL 23CM
HK 821-19	1	JANSEN-MIDDLETON SEPTUM FORCEPS 19CM
NC 300-14	2	CUSHING DURA HOOK 14CM
OD 942-21	1	HAJEK MALLET SOLID STEEL 200G 20CM